



“ADHESION WHEN YOU NEED IT.”

The reliable temporary eugenol-free cement with calcium hydroxide

- Safe adhesion
- Pulp protection
- Bacteriostatic

Provicol



Provicol

INTERNATIONALLY PROVEN. FOR YEARS.

Indications

Provicol is a eugenol-free temporary cement with calcium hydroxide for temporary luting of inlays, onlays, crowns, bridges, partial crowns, attachments, posts, and for temporary filling of small cavities.

Advantages

- Tough-elastic consistency for safe adhesion
- Calcium hydroxide
 - Supports the formation of secondary dentine
 - Bacteriostatic effect
- Eugenol-free: no allergy potential
- Ideal before adhesive luting with resin cements

Additional advantages of Provicol QM/ Provicol C:

- Fast application
- Bubble-free perfect mixing
- Direct application

Available in 3 versions:

- Provicol QM in the practical QuickMix syringe
- Provicol C in cartridges
- Provicol in tubes for hand mixing



Presentations

Provicol QM

REF 1074 pack 5 ml QuickMix syringe, mixing tips type 10

REF 2189 mixing tips type 10, 50 pieces

Provicol C

REF 1076 pack 2 x 65 g cartridges, mixing tips type 6

REF 2185 mixing tips type 6, 50 pieces

Provicol

REF 1075 pack 25 g base, 25 g catalyst

REF 1078 bulk pack 5 x 2 tubes, catalyst and base 25 g

VOCO GmbH
 P.O. Box 767
 27457 Cuxhaven
 Germany

Tel.: +49 (0) 4721-719-0
 Fax: +49 (0) 4721-719-140

info@voco.de
 www.voco.com

Available from:

